

## SPECIFICATIONS

**Model**

TNT1111

**Weight**

209 - 229 lbs

\*Intended for indoor use only

**Dimensions**

Length 86.1"

Width 38.5"

Height 37.1"

**Recommended Mattress**

36"x80"



## DETAILS

**Color**Steel Color: Khaki, EbonyOptional Steel Colors: Cream, London Fog, Pewter, Shark, Slate, Irish Moss, Dusk, Mink

\*Minimums apply

Laminate: Wild CherryOptional Laminate Colors: Kensington

Maple, Rustic Slate, Pewter Brush

\*Minimums apply

\*Due to differences in monitors, printers and materials, actual product colors may vary. Color chips are available for more accurate depiction of color.

**Options**

Steel Drawers Model TNT1012, Laminate Drawers Model TNT1212, Headboard &amp; footboard (Wild Cherry Standard), Ladder, Safety Rail, Bunk Lock Clamps, Glide Kit, Bolt-down Kit, Bottom Closure Kit

**Assembly**

Attach end frames to bed base with furnished hex socket button head fasteners (threaded inserts in vertical members of end frame). Optional accessories attach to bed (with hardware supplied through purchase of accessories).

**Fabrication & Material**

Cut, formed, welded steel, powder coated steel.

Bed deck: 18 gaugeBody: 18 gaugeGussets: 13 gaugeEnd frame legs: 16 gaugeEnd frame cross-members: 16 gaugeOptional Headboard: 45 lb density M3 grade particle board with durable high-pressure laminate surface and 2mm PVC edge (Wild Cherry Standard).

Made in the U.S.A.

**Installation**

None required for standard model

**Warranty**

10 year limited replacement warranty

**Maintenance**

Norix recommends the use of soap and water for both the regular cleaning maintenance and deep-scrub cleaning of powder coated steel.

Optional Decorative Laminate: Norix recommends the use of Simple Green® or an equivalent cleaning product.

## TESTING

**Environmental Test Standards**

- GREENGUARD and GREENGUARD Gold Certified by UL Environment. Product certified for low chemical emissions to UL 2818. UL.COM/GG.
- Compliant with Healthier Hospitals guidelines

**Performance Test Standards**

Tested to 800 lb static load